

beta Tubulin Monoclonal Antibody

Catalog Number: E-AB-20033 17 Publications



Note: Centrifuge before opening to ensure complete recovery of vial contents.


Description

Reactivity	Human, Mouse, Rat, Monkey, Chicken, Dog, Hamster, Rabbit, Sheep, Insect, Yeast
Immunogen	Synthetic Peptide
Host	Mouse
Isotype	IgG
Clone	Clone: 8B2
Purification	Protein A purification
Conjugation	Unconjugated
Formulation	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.

Applications Recommended Dilution


WB	1:5000-1:10000
IHC	1:100-1:300
IF	1:100-1:300

Data




Western Blot analysis of A549, Rat brain, Mouse brain, Chicken lung, Rabbit testis, Sheep muscle using beta Tubulin Monoclonal Antibody at dilution of 1:5000.

Observed Mw: 55kDa
Calculated Mw: 50kDa



Immunohistochemistry of paraffin-embedded Mouse testis tissue using beta Tubulin Monoclonal Antibody at dilution of 1:200.



Immunofluorescence analysis of HeLa tissue using beta Tubulin Monoclonal Antibody at dilution of 1:100.

Preparation & Storage

Storage Store at -20°C. Avoid freeze / thaw cycles.

For Research Use Only

A Reliable Research Partner in Life Science and Medicine

Toll-free: 1-888-852-8623

Web: www.elabscience.com

Tel: 1-832-243-6086

Email: techsupport@elabscience.com

Fax: 1-832-243-6017

beta Tubulin Monoclonal Antibody

Catalog Number: E-AB-20033 17 Publications



Background

There are five tubulins in human cells: alpha, beta, gamma, delta, and epsilon. Tubulins are conserved across species. They form heterodimers, which multimerize to form a microtubule filament. An alpha and beta tubulin heterodimer is the basic structural unit of microtubules. The heterodimer does not come apart, once formed. The alpha and beta tubulins, which are each about 55 kDa MW, are homologous but not identical. Alpha, beta, and gamma tubulins have all been used as loading controls. Tubulin expression may vary according to resistance to antimicrobial and antimitotic drugs.

For Research Use Only

A Reliable Research Partner in Life Science and Medicine

Toll-free: 1-888-852-8623

Web: www.elabscience.com

Tel: 1-832-243-6086

Email: techsupport@elabscience.com

Fax: 1-832-243-6017