

## Dulbecco's Phosphate Buffer (DPBS), without calcium, magnesium,

#### phenol red, 10×

Cat. No: PB180532 Size: 500mL

# General Information

Product Form	Liquid
Concentration	10×
pH Flabscienc	7.2-7.4
Concentration pH D-Glucose Elabscienc	negative
Phenol red	negative
CaCl2	negative
MgCl2•6H2O	negative
Sodium pyruvate	negative
Storage	2-30℃, Shading Light
Shipping	RT
Expiration date	24 months



## Background

Balanced Salt Solution (Physiological Solution) have the properties of buffer solution (regulate pH), normal saline (maintain osmotic pressure) and culture medium (provide nutrition). It can meet the basic needs of survival and metabolism of tissues, organs or cells in vitro. A small amount of phenolic red was added to some equilibrium salt solutions to indicate the pH change of the solution.

Dulbecco's Phosphate Buffered Saline (DPBS) is one of the most widely used equilibrium salt solutions in biochemistry. The main components are NaCl, KCl, KH2PO4 and Na2HPO4. DPBS can be divided into two types according to whether contain calcium and magnesium. Different from conventional PBS, DPBS phosphate content was slightly lower. DPBS was mainly used for embryological research.

# Notesby Elab

- 1. This product is for research use only.;
- 2. This product is sterilized by  $0.1 \mu m$  filtration.;
- 3. It is necessary to pay attention to the aseptic operation and avoid the contamination.