

## Elab Fluor® 594-Rabbit Anti-Human IgG(H+L)

catalog number: E-AB-1146

**Note:** Centrifuge before opening to ensure complete recovery of vial contents.

| Description           |                                                                                                             |
|-----------------------|-------------------------------------------------------------------------------------------------------------|
| Reactivity            | Human                                                                                                       |
| Immunogen             | Human IgG                                                                                                   |
| Host                  | Rabbit                                                                                                      |
| Purification          | Affinity purification                                                                                       |
| Conjugation           | Elab Fluor® 594(Ex/Em:591nm/614nm)                                                                          |
| Buffer                | PBS with 0.1% Proclin300, 1% protective protein and 50% glycerol, pH 7.4                                    |
| Applications          | Recommended Dilution                                                                                        |
| IF                    | 1:150-1:300                                                                                                 |
| FCM                   | 1:250-1:500                                                                                                 |
| Preparation & Storage |                                                                                                             |
| Storage               | Store at -20°C Valid for 12 months. Avoid freeze / thaw cycles. Protected from prolonged exposure to light. |
| Shipping              | Ice bag                                                                                                     |