The professional cell culture enpowers a healthier world

DMEM (Low glucose), with L-alanyl-L-glutamine, HEPES, without

sodium pyruvate

Cat. No: PM150263

Size: 500mL

General Information

Product Form Liquid Concentration 1×

pH 7.2-7.4

D-Glucose 1000mg/L

HEPES 25mM

L-Alanyl-L-Glutamine 4mM

NaHCO3 3700mg/L

Phenol red 15mg/L

Storage 2-8℃, Shading Light

Negative

Shipping RT

Expiration date 24 months



cella

Background

Sodium Pyruvate

DMEM (Dulbecco's Modified Eagle Medium) was developed on the basis of MEM medium. Compared with MEM medium, the content of amino acid increased by 2 times, the content of vitamin increased by 4 times, and the content of non-essential amino acid, trace iron ion and sodium pyruvate were increased by 4 times. The glucose content of DMEM medium was originally designed as 1000 mg/L (low Glucose type), and then developed into 4500mg/L (high Glucoser type), which has been widely used in cell culture.

DMEM (Low glucose) contains many kinds of amino acids, vitamins, inorganic salts and other ingredients for cell culture, but does not contain protein, lipids or any growth factors, so the product should be used with serum or serum-free additives.

Notes

- 1. This product is for research use only.;
- 2. This product is sterilized by 0.1µm filtration.;
- 3. It is necessary to pay attention to the aseptic operation and avoid the contamination.

Toll-free: 1-888-852-8623 Tel: 1-832-243-6086 Fax: 1-832-243-6017

Web: www.pri-cella.com Email: techsupport@pri-cella.com