

10× Normal goat Serum

Cat. No: E-IR-R111

Size: 5 mL / 10 mL

Cat.	Products	5 mL	10 mL	Storage
E-IR-R111	10× Normal goat Serum	5 mL	10 mL	-20°C

Introduction

This product is normal goat Blocking Buffer, added stabilizer, As a blocking buffer, It can effectively reduce the nonspecific staining of IHC/IF experiment, so as to reduce the background.

Instructions

The product is concentrated normal goat Blocking Buffer. Please dilute with PBS to 1× working solution before use.

For example: take 10 µL Normal Goat Blocking Buffer, add it to 90 µL PBS and mix.

Storage

Store at -20°C for 12 months. Avoid freeze / thaw cycles.

Cautions

1. This product couldn't be used for clinical diagnosis or treatment, food or medicine, and can't be stored in residence.
2. For your safety and health, please wear the lab coat and disposable gloves before the experiments.

For Research Use Only