

## Neuronal Medium (NM)

Cat. No.: CM-ZY003

Size: 500mL

### General Information

<b>Product Form</b>	Liquid
<b>Concentration</b>	1×
<b>Size</b>	500mL
<b>The proportion of complete medium preparation</b>	96% Basic Culture Medium for Neuron Cells + 2% B-27 Neuron Supplement, with vitamin A, 50× [PB180637] + 1% Supplement for Neuron Cells + 1% Additive for Neuron Cells
<b>Bacterial detection</b>	Negative
<b>Fungal detection</b>	Negative
<b>Mycoplasmal detection</b>	Negative
<b>Cells growth experiment</b>	The cells grow well and have a normal morphology
<b>Endotoxin level</b>	< 3 EU/mL
<b>Storage</b>	Store at the temperature indicated on the label, shading light
<b>Shipping</b>	Ice bag

### Background

NM has been meticulously optimized by the Procell's technical team. After long-term testing, this product can maintain the best growth state of neuron cells. This product already contains all the components required for the growth of neuron cells. If added in the recommended proportion, no other components need to be added and it can be directly used for the in vitro culture of neuron cells. This product is for further research use only and must not be used for diagnosis, treatment, clinical, household or any other purpose.

### Components of the culture medium

Components	Volume	Storage	Expiration date
Basic Culture Medium for Neuron Cells	500mL	2-8°C, shading light	12 months
B-27 Neuron Supplement, with vitamin A, 50×	10mL	-5~-20°C, shading light	12 months
Supplement for Neuron Cells	5mL	-5~-20°C, shading light	12 months
Additive for Neuron Cells	5mL	-5~-20°C, shading light	12 months

### Note

1. This product is only used for scientific research or further research, not for diagnosis and treatment.
2. It is necessary to pay attention to the aseptic operation and avoid the contamination during the culture.
3. B-27 Neuron Supplement, with vitamin A, 50×, Supplement for Neuron Cells and Additive for Neuron Cells should be stored at the temperature indicated on the label. It can be used after dissolving at 2-8°C. To avoid repeated freezing and thawing, it is recommended to store it in small quantities after subpackaging.
4. This product is packaged in a single component form. It needs to be prepared into a complete culture before use. The validity period of the complete culture medium is 3 months.
5. In order to maintain the best use effect of this product, it should not be placed in room temperature or high temperature environment for a long time.
6. All products should be used within their expiration dates. Once the expiration date is exceeded, they must be discarded.