



## DMEM (High glucose), ATCC Modification

Cat. No.: PM150293

Size: 500mL

## **General Information**

**Product Form** Liquid Concentration 1 × рH 7.2 - 7.4**D-Glucose** 4500 mg/L **HEPES** Negative 4 mM L-Glutamine Elabscience NaHCO<sub>3</sub> 1500 mg/L 15 mg/L Phenol red 1 mM Sodium pyruvate

2-8°C, Shading Light **Storage** ricella **Shipping** Room Temperature

**Expiration date** 24 months

## **Background**

Elabscience DMEM (Dulbecco's Modified Eagle Medium) was developed on the basis of MEM medium. Compared with MEM medium, the content of amino acid increased by 2 times, the content of vitamin increased by 4 times, and the content of non-essential amino acid, trace iron ion and sodium pyruvate were also increased. The glucose content of DMEM medium was originally designed as 1000 mg/L (low Glucose type), and then developed into 4500 mg/L (high Glucose type), which has been widely used in cell culture. DMEM (High glucose), ATCC Modification reduces the content of sodium bicarbonate.

This product contains many kinds of amino acids, vitamins, inorganic salts and other ingredients for cell culture, but does not contain protein, lipids or any growth factors, so the product should be used with serum or serum-free additives.

## Note

- This product is sterilized by 0.1 µm filtration, please pay attention to sterile operation when using it to avoid contamination;
- Do not freeze or thaw this product to maintain the best use effect of this product; 2.
- This product is only used for scientific research or further research, and is not used for diagnosis and PT1CEnot by Elabscience treatment.

Toll-free: 1-888-852-8623 Tel: 1-832-243-6086 Fax: 1-832-243-6017

Email: techsupport@pri-cella.com Web: www.pri-cella.com