icella®

Elabscience



## Ham's F-12, without L-glutamine

Cat. No.: PM150811

Size: 500mL

## **General Information**

**Product Form** Liquid Concentration 7.2-7.4 **D-Glucose** 1802 mg/L **Sodium Bicarbonate (NaHCO3)** 1176 mg/L L-Glutamine Negative **HEPES** Negative Phenol Red 1.2 mg/L **Sodium Pyruvate** 1 mM

2-8°C, Shading Light Storage

Shipping Room Temperature

**Expiration date** 24 months

## **Background**

Ham 's F-12 Nutrient Mixture was designed by Ham in 1969 based on Ham 's F-10 nutrient mixture, which was originally used for serum-free culture of CHO cells. Ham 's F-12 is often used as the basic culture medium for serum free culture. At low serum levels, it is especially suitable for single cell culture and clonal culture. After adding serum, it is also widely used in the culture of cancer cells and primary cells, such as rat hepatocytes, rat prostate epithelial cells, chondrocytes, rat myoblasts, and chicken embryonic cells. In addition, the DMEM/F12 medium obtained by mixing Ham 's F-12 in equal volume with DMEM is more nutritious and widely used. **Elabscience** 

Web: www.pri-cella.com

- 1. This product is for research use only.
- 2. This product is sterilized by 0.1 µm filtration.
- L'I Contam by Elabscience 3. It is necessary to pay attention to the aseptic operation and avoid the contamination during the culture.

Toll-free: 1-888-852-8623 Tel: 1-832-243-6086 Fax: 1-832-243-6017

Email: techsupport@pri-cella.com