

Recombinant Human SUMO3/SMT3A Protein (His Tag)

Catalog Number: PKSH033070

DIA • AN[®]
by Elabscience

Note: Centrifuge before opening to ensure complete recovery of vial contents.


Description

Synonyms	Small Ubiquitin-Related Modifier 3;SUMO-3;SMT3 Homolog 1;SUMO-2;Ubiquitin-Like Protein SMT3B;Smt3B;SUMO3;SMT3B;SMT3H1
Species	Human
Expression Host	E.coli
Sequence	Met 1-Phe103
Accession	P55854
Calculated Molecular Weight	13.8 kDa
Observed molecular weight	19 kDa
Tag	N-His

Properties

Purity	> 95 % as determined by reducing SDS-PAGE.
Endotoxin	< 1.0 EU per µg of the protein as determined by the LAL method.
Storage	Generally, lyophilized proteins are stable for up to 12 months when stored at -20 to -80°C. Reconstituted protein solution can be stored at 4-8°C for 2-7 days. Aliquots of reconstituted samples are stable at < -20°C for 3 months.
Shipping	This product is provided as lyophilized powder which is shipped with ice packs.
Formulation	Lyophilized from a 0.2 µm filtered solution of 20mM PB, 150mM NaCl, pH 7.2. Normally 5 % - 8 % trehalose, mannitol and 0.01% Tween80 are added as protectants before lyophilization. Please refer to the specific buffer information in the printed man
Reconstitution	Please refer to the printed manual for detailed information.

Data



> 95 % as determined by reducing SDS-PAGE.

Background

SUMO3 belongs to the SUMO protein family and operates like ubiquitin. Ubiquitin-like protein which can be covalently attached to target lysines either as a monomer or as a lysine-linked polymer. Nevertheless unlike ubiquitin that targets proteins for degradation; SUMO3 takes part in several cellular processes; such as nuclear transport; transcription regulation; apoptosis and protein stability. SUMO3 participates in amyloid beta generation and has a key role in the onset or progression of Alzheimer's disease.

For Research Use Only

A Reliable Research Partner in Life Science and Medicine

Toll-free: 1-888-852-8623

Web: www.elabscience.com

Tel: 1-832-243-6086

Email: techsupport@elabscience.com

Fax: 1-832-243-6017