

Recombinant Human SPINK1 Protein (His Tag)

Catalog Number:PKSH033046



Note: Centrifuge before opening to ensure complete recovery of vial contents.


Description

Synonyms	Pancreatic Secretory Trypsin Inhibitor;Serine Protease Inhibitor Kazal-Type 1;Tumor-Associated Trypsin Inhibitor;TATI;SPINK1;PSTI
Species	Human
Expression Host	HEK293 Cells
Sequence	Asp24-Cys79
Accession	P00995
Calculated Molecular Weight	7.3 kDa
Observed molecular weight	12-16 kDa
Tag	C-His

Properties

Purity	> 95 % as determined by reducing SDS-PAGE.
Endotoxin	< 1.0 EU per µg of the protein as determined by the LAL method.
Storage	Store at < -20°C, stable for 6 months. Please minimize freeze-thaw cycles.
Shipping	This product is provided as liquid. It is shipped at frozen temperature with blue ice/gel packs. Upon receipt, store it immediately at < - 20°C.
Formulation	Supplied as a 0.2 µm filtered solution of 20mM Tris-HCl, 500mM NaCl, 5% Trehalose, 5% Mannitol, 0.02% Tween 80, pH 9.0.
Reconstitution	Not Applicable

Data



> 95 % as determined by reducing SDS-PAGE.

Background

Serine Protease Inhibitor Kazal-Type 1 (SPINK1) is a trypsin inhibitor that prevent the trypsin-catalyzed premature activation of zymogens within the pancreas. Defects in SPINK1 are a cause of pancreatitis (PCTT). A disease characterized by the presence of calculi in pancreatic ducts. It causes severe abdominal pain attacks. Defects in SPINK1 are the cause of susceptibility to tropical calcific pancreatitis (TCP). Recombinant SPINK1 protein (rSPINK1) stimulated cell proliferation in benign RWPE as well as cancerous prostate cells. The research result indicated that the potential of SPINK1 as an extracellular therapeutic target in prostate cancer. In contrast, knockdown of SPINK1 in 22RV1 cells inhibited cell proliferation, cell invasion, and tumor growth in xenograft assays.

For Research Use Only

A Reliable Research Partner in Life Science and Medicine

Toll-free: 1-888-852-8623

Web: www.elabscience.com

Tel: 1-832-243-6086

Email: techsupport@elabscience.com

Fax: 1-832-243-6017