

Recombinant Human GADD45 β /GADD45B Protein (His Tag)

DIA • AN[®]
by Elabscience

Catalog Number:PKSH032467

Note: Centrifuge before opening to ensure complete recovery of vial contents.


Description

Synonyms	Growth Arrest and DNA Damage-Inducible Protein GADD45 Beta;Myeloid Differentiation Primary Response Protein MyD118;Negative Growth Regulatory Protein MyD118;GADD45B;MYD118
Species	Human
Expression Host	E.coli
Sequence	Met 1-Arg160
Accession	O75293
Calculated Molecular Weight	20.0 kDa
Observed molecular weight	24 kDa
Tag	N-His

Properties

Purity	> 95 % as determined by reducing SDS-PAGE.
Endotoxin	< 1.0 EU per μ g of the protein as determined by the LAL method.
Storage	Generally, lyophilized proteins are stable for up to 12 months when stored at -20 to -80°C. Reconstituted protein solution can be stored at 4-8°C for 2-7 days. Aliquots of reconstituted samples are stable at < -20°C for 3 months.
Shipping	This product is provided as lyophilized powder which is shipped with ice packs.
Formulation	Lyophilized from a 0.2 μ m filtered solution of 20mM Tris-HCl, 150mM NaCl, pH 8.0. Normally 5 % - 8 % trehalose, mannitol and 0.01% Tween80 are added as protectants before lyophilization. Please refer to the specific buffer information in the print
Reconstitution	Please refer to the printed manual for detailed information.

Data



> 95 % as determined by reducing SDS-PAGE.

Background

Growth Arrest and DNA Damage-Inducible Protein GADD45 β (GADD45B) is a member of the GADD45 family. GADD45B has been shown to interact with MAP3K4, ASK1, MAP2K7, and GADD45GIP1. GADD45B is involved in the regulation of growth and apoptosis. GADD45B reacts to environmental stresses by mediating activation of stress-responsive MTK1/MEKK4 kinase, which is an upstream activator of both p38 and JNK MAPKs. In addition, GADD45B participates in the down-regulation of hepatocellular carcinoma (HCC). It may serve as a possible therapeutic target.

For Research Use Only

A Reliable Research Partner in Life Science and Medicine

Toll-free: 1-888-852-8623

Web: www.elabscience.com

Tel: 1-832-243-6086

Email: techsupport@elabscience.com

Fax: 1-832-243-6017